



LXVB-01




THICKNESS 45MM

MODEL NUMBER & SIZE	LXVB-01 (SINGLE FILM)
PARAMETERS	
VIEW AREA (INCHES)	14 X 17"
EXTERNAL DIMENSIONS (MM)	496 X 512 X 47 (MM)
COLOR TEMPERATURE	> 12000K
LIGHT SOURCE	LED
LED LIFE	100,000 HOURS
POWER	AC 100-250V (50HZ)
OUTPUT	12V
POWER SUPPLY	40W

ADDITIONAL FEATURES: THE SAME PRODUCT IS AVAILABLE WITH SENSOR (AUTO ON/OFF) & ALSO WITHOUT SENSOR.




LXVB-02



THICKNESS 45MM

MODEL NUMBER & SIZE	LXVB-02 (TWO FILMS)
PARAMETERS	
VIEW AREA (INCHES)	28 X 17"
EXTERNAL DIMENSIONS (MM)	870 X 512 X 47 (MM)
COLOR TEMPERATURE	>12000K
LIGHT SOURCE	LED
LED LIFE	100,000 HOURS
POWER	AC 100-250V (50HZ)
OUTPUT	12V
POWER SUPPLY	80W

ADDITIONAL FEATURES: THE SAME PRODUCT IS AVAILABLE WITH SENSOR (AUTO ON/OFF) & ALSO WITHOUT SENSOR.




LXVB-03



THICKNESS 45MM

MODEL NUMBER & SIZE	LXVB-03 (THREE FILMS)
PARAMETERS	
VIEW AREA (INCHES)	42 X 17"
EXTERNAL DIMENSIONS (MM)	1217 X 512 X 47 (MM)
COLOR TEMPERATURE	>12000K
LIGHT SOURCE	LED
LED LIFE	100,000 HOURS
POWER	AC 100-250V (50HZ)
OUTPUT	12V
POWER SUPPLY	120W

ADDITIONAL FEATURES: THE SAME PRODUCT IS AVAILABLE WITH SENSOR (AUTO ON/OFF) & ALSO WITHOUT SENSOR



LXVB-04



THICKNESS 45MM

MODEL NUMBER & SIZE	LXVB-04 (FOUR FILMS)
PARAMETERS	
VIEW AREA (INCHES)	56 X 17"
EXTERNAL DIMENSIONS (MM)	1570 X 512 X 47 (MM)
COLOR TEMPERATURE	>12000K
LIGHT SOURCE	LED
LED LIFE	100,000 HOURS
POWER	AC 100-250V (50HZ)
OUTPUT	12V
POWER SUPPLY	160W

ADDITIONAL FEATURES: THE SAME PRODUCT IS AVAILABLE WITH SENSOR (AUTO ON/OFF) & ALSO WITHOUT SENSOR