

TONISHA
Lets Move Towards Visibility

*Opton 2004 represents high definition
Orthopantomograph Dental X-ray unit*

**BARC/AERB approved
ISO 9001 : 2000 COMPANY**

- * Standard Panoramic Projection
- * Child Projection
- * Tempo Mandibular Joints
- * Frontal Maxillary Sinuses
- * Cephalogram AP / Lateral View



AN ISO 9001 : 2000 COMPANY

Lateral Skull View




Panoramic View



OPTON 2004

TONISHA MEDICAL ORTHOPANORAMIC UNIT

Technical Data



Console



Sharp, Sleek Red Cross-Beam for Precise and accurate Positioning

Power Supply 230 V/50 Hz
Nominal Power 630 Watt
Max. Line Current 5 Amp.
Max. Power 1100 Watt
Displacement 100 cm.
Console/Generator Touch Panel, Two Line alphanumeric LCD Display, Two Pulse X-ray Generator

Remote Light weight, Easy operation for X-ray, Test Run and Reset.

Anode Voltage 55 Kvp to 85 Kvp (in 5 Kvp steps)
Anode Current 10 ma Fixed
X-ray tube CEI 105 (Made in ITALY)
Tube Focus .5mm x .5mm
Total Filtration 2.5mm Aluminium
Exposure time OPG: 17 Sec. Adult
 14 Sec. Child

Patient Adj. Light TMJ: 4x4 Sec.
Centering reference CEPH: 0.2 Sec. to 3 Sec.
Cassette Type Red laser Cross light
Cassette Size Chin rest with Bite wing
Vertical Movement Flat (General / Digital)
Cephalostat 6"x12", 8"x10", 10"x12
 Motorised
 AP/Lateral/Oblique View

- * Interlocking for all required Consents
- * Choice of Adult/Child programs (Three patient sizes)
- * The unit can be upgraded to Digital Model in future at site.

TONISHA
Lets Move Towards Visibility

TONISHA ELECTRONICS CORPORATION

Head Office : E-100, 1st Floor, Naraina,
New Delhi - 110028, (INDIA).
Works : B-222 (Rear Side), Naraina Industrial Area,
Phase-1, New Delhi - 110028, (INDIA).



SALES & SERVICE

UPASANA ENTERPRISES

B-303, Patidar Complex, Near Police Station,
Kannamwar Nagar - 2, Vikhroli (E), Mumbai - 83.
Tel.: 2577 4842 + Telefax: 2577 4879
Mobile: 98211 21823 / 98694 73756 / 98215 26497
E-mail: upasanaenterprises@yahoo.com