

TONISHA

Lets move towards visibility

High definition Panoramic Dental X-Ray unit from the Professional Company

BARC/AERB approved

ISO 9001 : 2000 COMPANY



- *Standard Panoramic Projection*
- *Child Projection*
- *Tempo Mandibular Joints*
- *Frontal Maxillary Sinuses*
- *Cephalogram AP / Lateral View*

Lateral Skull View



Panoramic View



OPTON 2004

Flat Cassette Orthopantomograph unit designed for optimal Patient positioning & High quality Images

TONISHA MEDICAL ORTHOPANORAMIC UNIT

Technical Data



Console



Sharp, Sleek Red Cross-Beam for Precise and accurate Positioning

Power Supply	230 V/50 Hz
Nominal Power	630 Watt
Max. Line Current	5 Amp.
Max. Power	1100 Watt
Displacement	100 cm.
Console/Generator	Touch Panel, Two Line alphanumeric LCD Display, Two Pulse X-ray Generator

Remote Light weight, Easy operation for X-ray, Test Run and Reset.

Anode Voltage 55 Kvp to 85 Kvp (in 5 Kvp steps)

Anode Current 10 ma Fixed
X-ray tube CEI 105 (Made in ITALY)

Tube Focus .5mm x .5mm

Total Filtration 2.5mm Aluminium

Exposure time OPG: 17 Sec. Adult
14 Sec. Child

TMJ: 4x4 Sec.
CEPH: 0.2 Sec. to 3 Sec.

Patient Adj. Light Red laser Cross light
Centering reference Chin rest with Bite wing

Cassette Type Flat (General / Digital)

Cassette Size 6"x12", 8"x10", 10"x12

Vertical Movement Motorised

Cephalostat AP/Lateral/Oblique View

**Interlocking for all required Consents*

**Choice of Adult/Child programs (Three patient sizes)*

**The unit can be upgraded to Digital Model in future at site.*

TONISHA

Lets Move Towards Visibility

TONISHA ELECTRONICS CORPORATION

Head Office : E-100, 1st Floor, Naraina,

New Delhi - 110028, (INDIA).

Works : B-222 (Rear Side), Naraina Industrial Area,

Phase-1, New Delhi - 110028, (INDIA),



SALES & SERVICE

UPASANA ENTERPRISES

B-303, Patidar Complex, Near Police Station,
Kannamwar Nagar - 2, Vikhroli (E), Mumbai - 83.

Tel.: 2577 4842 • Telefax: 2577 4879

Mobile: 98211 21823 / 98694 73756 / 98215 26497

E-mail: upasanaenterprises@yahoo.com