

NEW

Reversible Apron

Techset Reversible Apron provides two different colours to facilitate wearing both sides according to your mood and environment.

Available with / without attached thyroid collar.



Coat Apron (Velcro Type)

Techset Coat Apron is stylish and provides the highest possible levels of comfort without compromising on quality.

Its comfortable design prevents fatigue for prolonged periods of time.



Economy Apron (Belt Type)

Techset Economy Apron is economically designed to provide comfort without compromising on quality.

Its belt-type fitting provides a unique flexibility and durability to the product.


Wrap Around Apron (Double Sided)

Convenience is another name for the **Techset Wrap Around**. Apart from providing all-round protection, it also fits the body comfortably enough to help prevent discomfort when worn for longer periods of time.



Skirt & Vest


Techset Skirt-Vest duo is a 2-piece protective gear for those who wish to be clad in style. It is a hot favourite among radiologists and their colleagues who have to keep moving in and out of the X-ray room.



Ovarian Shield



Gonad Shield



Techset Ovarian & Gonad Shield ensures that the male genital region and female ovarian region is adequately shielded against X-ray penetration.

Care has been taken to avoid irritation to men wearing Gonad Shield & women wearing Ovarian Shield.

Mini Apron

When you require protection only below the abdomen, **Techset Mini Apron** scores a notch above others.

Stylish, elegant and light-weight are of the attributes of the Mini Apron that make it a preferred choice of radiologists.



Dental Apron

Techset Dental Apron with attached thyroid collar is the best X-ray protection that a dentist can give to patients.

Light-weight & comfortable, assures patient's safety while helping you bid good-bye to your radiation worries.

Thyroid Shield

Techset Thyroid Shield is the ultimate protection for the neck against X-ray radiation. While its neck-high design ensures complete covering upto the chin, its Velcro-belt offers the flexibility to fit any neck size.



Our strengths lie in our passion for research and the ability to bring technology to the people at affordable rates. No wonder, Techset is today a sought-after name among the prominent radiologists.

Our speciality is comfortable, safe, light-weight, flexible and extremely crack-resistant Lead Vinyl Aprons providing the best available protection against X-ray radiations.

Our products can be customized as per your requirement within your time frame.

Standard Sizes

Coat Apron and Wrap Around Apron

Sizes	Length	Width
Regular	100 cm	60 cm
Medium	110 cm	60 cm
Long	120 cm	60 cm
Short	90 cm	60 cm

Skirt & Vest

	Sizes	Length	Girth
Skirt	Regular	60 cm	120 cm
	Long	70 cm	120 cm
Vest	Regular	70 cm	120 cm
	Long	80 cm	120 cm

Lead Equivalence

All our products are provided in 0.5 mm Pb.
Other lead equivalence available on request.

- ✓ **Extra Smooth Surface**
Reduced damage & longer life due to reduced friction
- ✓ **Adjustable Velcro Belts**
Comfortable fitting & easy adjustability to persons of any built & size
- ✓ **Easily Washable Material**
Makes the aprons / shields reusable while reducing the risk of cross infections
- ✓ **Evenly Distributed Lead Compound**
High level of comfort, uniform protection, flexibility and lightness
- ✓ **Multilayer Lead Vinyl Material**
Enhanced protection from X-ray matching AERB specifications radiations
- ✓ **Quality Check at Every Step**
Ensures durability and hassle-free use, making it one of the best products in the market

Techset has been a global player exclusively focusing on manufacturing of X-Ray Protection equipment since 1990. Conforming to international quality standards & manufacturing processes, and years of experience has made Techset a name you can trust.

Contact For Sales & Service **UPASANA ENTERPRISES**

Mob.: 98211 21823/ 98694 73756 / 98215 26497

Tel.: 022-2574 0383 / 2574 0112

Email.: upasanaenterprises@yahoo.com

Web.: www.upasanaenterprises.com

THE **ONLY** LEAD VINYL APPAREL Mfgd. in India for



X RAY PROTECTION

Techset



ISO 9001:2000
06L-0041